


Gelenkbolzenschelle



Größentabelle und Technische Daten

Artikel-Nr.	Größe Zoll	Spannbereich (mm)	VPE Box	Gewicht per box (kg)	b (mm)	Schraube	Schlüsselweite	s (mm)
TEK-01		13-15	50	2,2	18	M6	11	0,8
TEK-02		15-17	50	2,2				
TEK-03		17-19	50	2,2				
TEK-04		19-21	50	2,4				
TEK-05		21-23	50	2,5				
TEK-06		23-25	50	2,5				
TEK-07		25-27	50	2,6				
TEK-08		27-29	50	2,6				
TEK-09		29-31	50	2,8				
TEK-10	1"	31-33	50	2,8				
TEK-11		34-36	50	5,6	25	M8	14	1,0
TEK-12	1 1/4"	37-40	50	5,8				
TEK-13		41-43	50	6				
TEK-14	1 1/2"	44-46	50	6				
TEK-15		47-51	50	6,2				
TEK-16		52-55	50	6,4				
TEK-17	2"	56-61	50	6,4				
TEK-18		59-65	50	6,5				
TEK-19		62-67	50	6,6				
TEK-20		65-71	50	6,8				
TEK-21	2 1/2"	68-74	50	7				
TEK-22		72-78	50	7,2				
TEK-23		75-83	50	7,4				
TEK-24	3"	84-91	50	7,6				
TEK-25	3 1/2"	92-102	25	4				
TEK-26		100-108	25	4,4				
TEK-27		103-110	25	4,6				
TEK-28	4"	110-120	25	4,6				
TEK-29		121-131	25	4,6				
TEK-30	5"	135-145	25	4,8				
TEK-31		145-155	25	4,9				
TEK-32		155-165	25	5				
TEK-33	6"	165-175	25	5,1				



Material

W1	Band & Gehäuse	: 1.0330 [DC01 (C590) Qualitäts Stahl /
	Keine Verwendung von Chrom (VI)]	
	Schraube	: 1.5525 [(20MnB4) Spezialbeschichtet]
	Salzwassersprüh test:	200 Stunden
W2	Band & Gehäuse	: 1.4016 [(X6Cr17) Magneticst/st]
	Schraube	: 1.5525 [(20MnB4) Spezialbeschichtet]
	Salzwassersprüh test:	400 Stunden

